



TETRA**GRAPH**

ADULT
Sensor Placement

Quantitative
Train-of-Four Monitor



Ulnar Nerve

Option 1



Grounding

Landmark

Hand

1 cm Below
the Wrist Crease

Foot

Medial Malleolus

Option 2



AP*

Posterior Tibial Nerve

Option 3



FHB*

Stimulating Ulnar Nerve | Posterior Tibial Nerve



Visit [Senzime.com/Support](https://www.senzime.com/support)

For more information, scan the QR code.

*ADM Abductor Digiti Minimi muscle | AP Adductor Pollicis muscle | FHB Flexor Hallucis Brevis muscle

