




# Fo<sub>2</sub>cus Flo<sub>2</sub>w Open Oxygen Masks

## Fo<sub>2</sub>cus Flo<sub>2</sub>w



**Mask size now engraved on each mask!**  
**Color coding on the label for easy size identification!**

*Fo<sub>2</sub>cus Flo<sub>2</sub>w design concentrates the oxygen to the upper lip and nasal passages to ensure the maximum oxygen is delivered to the patients inspiratory pathway directly without dispersion into the full mask area as with normal closed oxygen masks. Positioning of mask is important to obtain the best FiO<sub>2</sub> levels %, 6 choices of mask are offered to allow best fitting options to patients. This also aids the performance of the mask with more size options.*

Adult sizes	Flow Rate	FiO <sub>2</sub> %
	1 LPM	24 - 26
	2 LPM	26 - 30
	3 LPM	30 - 35
	4 LPM	35 - 40
	5 LPM	38 - 50
	8 LPM	40 - 66
	10 LPM	50 - 75
	12 LPM	55 - 75
	15 LPM FLUSH	75 - 90
Child / Infant Sizes	Flow Rate	FiO <sub>2</sub> %
	0.25 LPM	22
	0.50 LPM	24
	0.75 LPM	27
	1 LPM	35
	1.5 LPM	40
	2 LPM	45
	3 LPM	48
	5+ LPM	60 - 65

**\*\* The Actual FiO<sub>2</sub> will vary with Patient's Rate & Depth of breathing, and the position / distance of patient's face to the oxygen inlet flow. As this affects the target area for the oxygen flow. Pulse oximeter use is recommended when deciding the correct flow rate for the patient.**

**With "ELLIE" Masks careful increase of the flow rate from 0.25 LPM upwards is recommended to achieve the desired flow rate needed. Again this is in conjunction with a Pulse Oximeter to verify actual FiO<sub>2</sub>**

### KEY FEATURES

- 100% NON PVC ECO Product - 6 Sizes: Adult XL, Large, Medium, Small, Child, Infant
- More sizes allow more precise fit and more accurate oxygen delivery
- Open design less claustrophobic & allows speech to be heard clearly
- Easy swivel oxygen tube connector allows movement without dislodging mask on patient face

REF CODE	DESCRIPTION	Pack Count Per Carton
22OP1441	FF X-Large Adult Mask with 2.1m Tubing	25
22OP1431	FF Large Adult Mask with 2.1m Tubing	25
22OP1421	FF Medium Adult Mask with 2.1m Tubing	25
22OP1411	FF Small Adult Mask with 2.1m Tubing	25
22OP1471	FF Child (ELLIE) Mask with 2.1m Tubing	25
22OP1481	FF Infant (ELLIE) Mask with 2.1m Tubing	25
22OP142130	FF Medium Adult Mask with 3.0m Tubing	25

# Product Technical Summary

100% PVC FREE

## Regulatory

[CONTACT US](#)

The product complies with legal requirements for medical devices ISO 11135



**REACH regulation:** The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

<b>Biological evaluation:</b>	EN ISO 10993-5, EN ISO 10993-7 and BS EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Devices Directive 93/42/EEC, Class IIa
<b>Compliance with product standards</b>	EN ISO 11135, EN ISO 11737-1, BS EN ISO 14971
<b>Quality standards</b>	ISO 13485
<b>Label information and packaging compliance</b>	EN ISO 15223-1, EN 1041 and EN ISO 11607-1 , BS EN 15986

## Materials: Mask TPO, Tubing TPE, Swivel Connector PP, Strap TPR

### Warnings and Recommendations

<b>Warnings</b>	The product is single use. If reused performance of product may deteriorate, cross contamination may occur. Single-patient-use longterm approved with cleanliness monitored.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	4 Years
<b>Country of origin</b>	UK
<b>Product Weight</b>	15 - 20g

### Packaging Information

Individual Pack	Material of Pack	Size (L x W x H) mm	Shipping Box	Size (L x W x H) mm
Single Pack Pouch	PE Bag	100 x 200 mm	Cardboard	L x W x H cm

# Barcodes / GTINs

REF	Single Pack	Transport Carton
22OP1441	<p>22OP1441 SINGLE:  (01) 0 5055276 74432 0</p>	 (01) 1 5055276 74432 7
22OP1431	<p>22OP1431 SINGLE:  (01) 0 5055276 74431 3</p>	 (01) 1 5055276 74431 0
22OP1421	<p>22OP1421 SINGLE:  (01) 0 5055276 74430 6</p>	<p>BOX:  (01) 1 5055276 74430 3</p>
22OP1411	 (01) 0 5055276 74471 9	 (01) 1 5055276 74471 6
22OP1471	<p>22OP1471 SINGLE:  (01) 0 5055276 74433 7</p>	 (01) 1 5055276 74433 4
22OP1481	<p>22OP1481 SINGLE:  (01) 0 5055276 74434 4</p>	 (01) 1 5055276 74434 1

Quality | Value | Integrity

Meditech Systems Limited

Shrublands Estate, Shearstock, Shaftesbury Dorset, SP7 9PT

[www.meditechsystems.co.uk](http://www.meditechsystems.co.uk)

Tel: +44 (0) 1747 821546 Email: [sales@meditechsystems.co.uk](mailto:sales@meditechsystems.co.uk)

facebook: [@meditechsystemslimited](https://www.facebook.com/meditechsystemslimited) twitter: [@MeditechLtd](https://twitter.com/MeditechLtd) LinkedIn: [meditech-systems-limited](https://www.linkedin.com/company/meditech-systems-limited)

© 2020 Meditech Systems Ltd | All rights reserved | FDA registration No.10048739