



4-year guarantee

# PARI SINUS2

## Technical datasheet

Precise, effective and gentle treatment of sinusitis.



**Technical data:**

**Alternating current:** 220 – 240V, 50Hz, **Pressure:** 1.5 bar, **Weight:** 1.7 kg, **Housing dimensions:** 18.5x13x15 cm (WxHxD), **Fill volume:** min. 2 ml, max. 8 ml, **Noise volume:** 55 db(A)

**Second to none** – pulsating aerosol effectively deposits the medicine or saline solution in the paranasal sinuses <sup>1, 2</sup>



Gamma camera images overlap with MRI images

**The PARI SINUS2 includes:**

- PARI LC SPRINT SINUS nebuliser and PARI SINUS Tube System



	NASAL SPRAY	PARI SINUS2
<b>Nasal deposition in %</b>	96.5*	21.5*
<b>Deposition in the paranasal cavities in %</b>	0.45*	12.5*

◀ **Around 28 times more effective than the nasal spray**

\* Deposition (in %) of the applied active agent.

1. Möller W. et al (2009): Ventilation and aerosolized drug delivery to the paranasal sinuses using pulsating airflow – a preliminary study. *Rhinology*, 47, 405-412, 2009.  
 2. Schuschnig U. et al (2006): Comparison of delivery efficiency in a nasal case model of fluticasone propionate suspensions and a novel solution aerosolized via the PARI Vibrent™.

# PARI SINUS2

## Item numbers for your region

Product	Sales Region	IFU languages	REF (Item No.)
PARI SINUS2 Nebuliser System 	Western Europe	de, fr, it, pt, es, nl, en	128G1001
	Eastern Europe	cs, pl, da, fi, no, sv	128G1006
	UK & Ireland	en	128G1004
	Australia & New Zealand	en	128G1015
PARI SINUS2 Year Pack 	Western Europe	de, fr, it, pt, es, nl, en	023G3811
	Eastern Europe	cs, el, hu, pl, da, fi, no, sv	023G3816

### Useful accessories for PARI SINUS2

Product	REF (Item No.)
PARI LC SPRINT SINUS Application Set	023G2840
PARI SINUS Tube System	041G4570
PARI LC SINUS Nose Plug (3 pieces)	041G0576
PARI LC SPRINT Nebuliser (WE) (IFU languages: de, fr, it, pt, es, nl, en)	023G1001
PARI LC SPRINT Nebuliser (EE) (IFU languages: bg, cs, el, et, hu, pl, ro, sk, sl)	023G1006
PARI Bag	130G5000



PARI GmbH  
 Moosstraße 3, 82319 Starnberg, GERMANY  
 E-Mail: [export@pari.de](mailto:export@pari.de)  
 Phone: +49 (0) 8151 279 220  
[www.pari.com](http://www.pari.com)